

Specifications

DPM2 PLUS Universal Pressure Meter

Pressure Measurement

DPM2 PLUS:

Operating range: -698 to +802 mmHg
-949 to +1090 cmH₂O
-374 to +429 inH₂O
-13.50 to +15.50 PSI
-13.5 to +100 PSI

Accuracy: ±1% of range

Units of measure: mmHg, cmH₂O, PSI, inH₂O

General Information

Temperature Requirements:	Operating: 0°C to +30°C (+32°F to +86°F)
Display/Control:	0.5 in LCD with LO BATT and negative polarity indication
Data Outputs:	0.1 V/psi (all ranges except 100 psi) 0.01 V/psi for 100 psi range
Power:	9 Volt Alkaline Battery
Dimensions:	H x W x D: 14.61cm x 9.14cm x 4.83cm (5.75 in x 3.6 in x 1.9 in)

Fluke Biomedical

PO Box 9090, Everett, WA 98206-9090 USA

Toll Free: (800) 648-7952, Fax (425) 446-5629, Email: sales@flukebiomedical.com

www.flukebiomedical.com