

Product Comparison

Pressure Meter Comparison Chart

Model	DPM3	DPM4-1H	DPM4-1G (replaces DPM-3)	DPM4-2H	DPM4-2G
Part number	Discontinued	2583121	2631330	2637760	2637772
Price	\$1,295	\$1,295	\$1,295	\$1,595	\$1,595
Operating range	-698 to +776 mmHg	-350 to +350 mmHg	-700 to +5000 mmHg	-350 to +350 mmHg	-700 to +5000 mmHg
Temperature simulation	Yes	No	No	No	No
Units of measure	cmH ₂ O, inH ₂ O, kPa, psi, mmHg	mBar, cmH ₂ O, psi, inHg, inH ₂ O, kgcm ² kPa, mmHg	mBar, cmH ₂ O, psi, inHg, inH ₂ O, kgcm ² kPa, mmHg	mBar, cmH ₂ O, psi, inHg, inH ₂ O, kgcm ² kPa, mmHg	mBar, cmH ₂ O, psi, inHg, inH ₂ O, kgcm ² kPa, mmHg
Temperature range	0 to +50 °C	-200 to +750 °C	-200 to +750 °C	-200 to +750 °C	-200 to +750 °C
Gas flow	None	None	None	-750 to +750 ml/min	-750 to +750 ml/min
Relative humidity	None	None	None	0 to 100% RH	0 to 100% RH
Barometric pressure	None	None	None	0 to 100% RH	0 to 100% RH

Fluke Biomedical

PO Box 9090, Everett, WA 98206-9090 USA

Tel: 425.446.6945, Fax: 425.446.5629, Email: sales@flukebiomedical.com

www.flukebiomedical.com