



Intermittent Suction Unit (ISU)

- ◆ Begins in the “OFF” Mode
- ◆ Six turns to full wall vacuum
- ◆ Front mounted on/off switch
- ◆ Glow-in-the-dark gauge increments and needle
- ◆ Source port at bottom of unit



Patented Push-To-Set™ Built-in feature eliminates the need to occlude patient tubing to set accurate vacuum levels

Push-To-Set™ Intermittent Suction Unit (PTS-ISU)

- ◆ Stat Mode starts in the “ON” Mode
- ◆ Quick-To-Max, Only 2 turns to full wall vacuum
- ◆ Top mounted on/off switch for ease of use
- ◆ Brighter glow-in-the-dark gauge increments and needle
- ◆ Source port at center of unit for space constraints
- ◆ Smaller, sleek design for tight spaces 6.5”H x 2.8”W x 4.8”D