

MAQUET

**OPERATING TABLE ACCESSORIES
NEW PRODUCTS AND INNOVATIONS**

SURGICAL WORKPLACES



ALWAYS STATE OF THE ART MAQUET ACCESSORIES

Tried and tested under practical conditions and daily use. MAQUET develops functional operating table accessories for optimum patient positioning. Ongoing development ensures that the high quality requirements can also be fulfilled in the future. The following excerpt from the complete MAQUET

accessories programme introduces the accessories that are in use every day. For a more comprehensive overview, you will find some of the previous models that are no longer available located next to the current products.

HEAD



1002.65A0

Connection fixture with central square socket

- Adjustable at three joints
- For mounting horseshoe-shaped head rests/plates and skull clamp
- Optimum adjustment range and safe positioning



1002.64A0

No longer available



1002.71A0

Horseshoe-shaped head rest, one-piece with fastening strap

- Detachable SFC pad
- Ideal pressure distribution in the supine and prone position
- For example, for connection fixture 1002.65A0



1002.6000

No longer available



1002.72A0

Horseshoe-shaped head rest, two-piece

- Head cup, width-adjustable
- Detachable SFC pad
- Ideal pressure distribution in the supine and prone position
- For example, for connection fixture 1002.65A0



1002.7000

No longer available



SAFE, RELIABLE, EFFICIENT: ACCESSORIES FOR THE DAILY CLINICAL ROUTINE MAQUET – THE GOLD STANDARD



OR equipment designed to meet the highest of demands.

MAQUET is one of the world leaders in advanced medical technology in OR and intensive care. Be it surgical workplaces, critical care or cardiovascular surgery. The needs of both personnel and patients take top priority when it comes developing products. MAQUET provides optimised process flows to meet the increasing demands made of safety, efficiency and cost effectiveness in the daily clinical routine.

In daily use: The appropriate accessories. A comprehensive range of operating table accessories completes the MAQUET portfolio. A selection of the most important currently available accessories can be found on the following pages. Of course, all products are also included in the complete accessories catalogue.

MAQUET – The Gold Standard.

**1002.73A0****Head plate, flat**

- Detachable SFC pad
- Support frame with integrated clamp, electrically conductive
- Optimised neck support
- Radiotranslucent in a.p. view in conjunction with 1005.5800

**1002.6100****No longer available****1002.05A0****Head support for shoulder surgery in beach chair position**

- Detachable SFC pad
- Two Velcro straps
- Safe positioning and lateral fixation
- For head plate adapter of shoulder plate
- Adjustable to head size
- Helmet, mechanical decontamination

**1001.44C0/D0****Arm posturing device**

- Side rail with fastening clamp, can be horizontally swivelled
- Arm board can be rotated, longitudinally and height-adjustable
- Arm rest (L 300 mm or 450 mm) adjustable

**1001.44.E0/F0****Arm posturing device**

- Like 1001.44C0/D0, but with detachable SFC pad
- Excellent pressure distribution

**1002.25A0****Arm protector, stainless steel**


- Detachable pad (L 460 mm, H 100 mm)
- Safe arm and patient positioning
- Integrated clamp for side rail attachment
- With 2 Velcro straps

**1002.17A0****Lateral support**


- Two PUR pads
- For pelvis support in lateral position
- Fastening to side rail with integrated clamp

**1001.56A0****Body strap with retaining clip**

- For fastening to side rails
- Washable
- Also for over-weight patients (W 120 mm, L 1500 mm)

**1001.56B0****Body strap with Velcro**

- For fastening to side rails
- Washable
- Also for over-weight patients (W 120 mm, L 1500 mm)

**1150.67BC****Pair of leg plates**

- Carbon-fibre version for 360° radioscopies
- Detachable SFC pad
- Can be abducted
- Length: 910 mm

1133.67BC**Pair of leg plates**

- As 1150.67BC, but matching mounting point 1130

**1001.65A0****Goepel knee crutch**

- 1 piece
- Knee crutch adjustable via stabilised ball-and-socket joint
- With pad and strap
- Adjustment via radial setting clamp (e. g. 1003.23C0)

**1001.73A0****Knee crutch, extra wide**

- For over-weight patients
- Detachable SFC pad
- Electrically conductive cover
- Fastening via radial setting clamp, e. g. 1003.23C0

**1002.55A0****Tube holder**

- For secure fixing of resuscitation tubes and monitor cables
- Rod with flexible shaft, height-adjustable
- Clamp must be ordered separately (e. g. 1003.23C0 or 1003.61A0)

**1009.01C0****Infusion stand**

- For fastening to the side rail using the integrated clamp
- Vertical stand, height-adjustable from 800-1500 mm or continuously adjustable via self-locking ball-and-socket joint

**1003.23C0****Radial setting clamp, can be screwed into place**

- For fastening various OR accessories to the side rail
- For rod diameters from 16 to 18 mm
- Locking of accessories and radial adjustment via tommy screw

**1003.61A0****Radial setting clamp, adaptable using the quick-action connector**

- Ideal clamping mechanism
- Radial adjustment through 360° for rod diameters of 16 to 18 mm

For local contact:

Please visit our Website
www.maquet.com

MAQUET

MAQUET GmbH & Co. KG
Kehler Straße 31
D-76437 Rastatt, Germany
Phone: +49 (0) 7222 932-0
Fax: +49 (0) 7222 932-571
info.sales@maquet.de
www.maquet.com

GETINGE

GETINGE Group is a leading global provider of equipment and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient handling and hygiene, disinfection, DVT prevention, medical beds, therapeutic surfaces and diagnostics. GETINGE provides solutions for infection control and prevention within healthcare and life sciences. MAQUET specializes in therapeutic applications, products, solutions and services for OR and ICU.