

Enteralite Infinity Administration Sets

Administration Sets

FOR THE ENTERALITE INFINITY ENTERAL FEEDING PUMP

Moog's line of Zevex Enteral feeding products include unique delivery sets that provide accurate delivery for the patient and ease-of-use for the caregiver.

Our EnteraLite Infinity enteral delivery sets come with pre-attached pump set, patented in-line occluder for automatic free-flow protection, built-in hangers, and a large top-fill opening with leak-proof cap to minimize formula spills and waste. All Moog enteral delivery sets are DEHP and latex-free.

INF1200


ENTERALITE INFINITY 1200 ML DELIVERY SET

The EnteraLite Infinity 1200 ml delivery set (INF1200) comes with pre-attached pump set, patented in-line occluder for automatic free-flow protection, built-in hanger, and a large top-fill opening with leak-proof cap to minimize formula spills and waste.

All Moog enteral delivery sets are DEHP and latex-free.

Non-IV Compatible

All mating adapters meet AAMI guidelines to prevent potentially hazardous situations in connecting to the patient's feeding tube.



Specifications	
Order Number	INF1200
Materials	DEHP and latex-free
Free-flow Protection	Patented in-line occluder for automatic free-flow protection

INF0500


ENTERALITE INFINITY 500 ML DELIVERY SET

The EnteraLite Infinity 500 ml delivery set (INF0500) comes with pre-attached pump set, patented in-line occluder for automatic free-flow protection, built-in hanger, and a large top-fill opening with leak-proof cap to minimize formula spills and waste.

All Moog enteral delivery sets are DEHP and latex-free.

Non-IV Compatible

All mating adapters meet AAMI guidelines to prevent potentially hazardous situations in connecting to the patient's feeding tube.



Specifications	
Order Number	INF0500
Materials	DEHP and latex-free
Free-flow Protection	Patented in-line occluder for automatic free-flow protection

INF0010

ENTERALITE INFINITY SPIKE SET

The EnteraLite Infinity spike set (INF0010) is used with pre-filled containers to reduce nursing time and minimize the risk of bacterial contamination. The INF0010 is non-IV compatible, includes a "for enteral feeding only" tag, comes with pre-attached pump set, and has a patented in-line occluder for automatic free-flow protection.

All Moog enteral delivery sets are DEHP and latex-free.

Non-IV Compatible

All mating adapters meet AAMI guidelines to prevent potentially hazardous situations in connecting to the patient's feeding tube.



Specifications	
Order Number	INF0010
Materials	DEHP and latex-free
Free-flow Protection	Patented in-line occluder for automatic free-flow protection