

## Enteral Accessories

### Mini Backpack

FOR BOTH ENTERALITE INFINITY AND ENTERALITE ENTERAL FEEDING PUMPS



This versatile pack can hold either the EnteraLite Infinity or EnteraLite enteral feeding pump with a 500ml bag in the front section and can also accommodate a 1200ml bag in the rear section. Also included in this pack is a pocket that may be used to hold an ice pack.

Specifications	
Order Numbers	PCK1001 - Red, PCK1002 - Green/Gray, PCK1003 - Black
Feeding Pump Compatibility	Works with either EnteraLite Infinity or EnteraLite enteral feeding pumps
Delivery Set Compatibility	Can accommodate either 500 ml or 1200 ml delivery sets
Dimensions	13" (33 cm) H x 8" (20.3 cm) W x 4.5" (11.4 cm) D

### Super-Mini Backpack

FOR BOTH ENTERALITE INFINITY AND ENTERALITE ENTERAL FEEDING PUMPS



Ideally suited for young pediatric patients, this tiny pack holds either the EnteraLite Infinity or EnteraLite enteral feeding pump with a 500ml bag.

Specifications	
Order Numbers	PCK2001 - Black, PCK2002 - Green/Gray, PCK2003 - Specialty
Feeding Pump Compatibility	Works with either EnteraLite Infinity or EnteraLite enteral feeding pumps
Delivery Set Compatibility	500 ml delivery set only
Dimensions	9.5" (24.1 cm) H x 8" (20.3 cm) W x 4" (10.2 cm) D

## Enteral Accessories

### Waist Pack

FOR BOTH ENTERALITE INFINITY AND ENTERALITE ENTERAL FEEDING PUMPS



Well-suited for adolescents and adults, this small pack holds either the EnteraLite Infinity or EnteraLite enteral feeding pump with a 500ml bag and can be worn on a belt at the hip, as a waist pack, or as a shoulder bag.

Specifications	
Order Numbers	PCK3001
Feeding Pump Compatibility	Works with either EnteraLite Infinity or EnteraLite enteral feeding pumps
Delivery Set Compatibility	500 ml delivery set only
Dimensions	7" (17.8 cm) H x 11" (27.9 cm) W x 3" (7.6 cm) D

### Backpack

FOR BOTH ENTERALITE INFINITY AND ENTERALITE ENTERAL FEEDING PUMPS




This larger, adult-sized pack can hold either the EnteraLite Infinity or EnteraLite enteral feeding pump with either a 500 ml or 1200 ml feeding bag. Also included in this pack is a pocket that may be used to hold an ice pack.

Specifications	
Order Numbers	PCK4001
Feeding Pump Compatibility	Works with either EnteraLite Infinity or EnteraLite enteral feeding pumps
Delivery Set Compatibility	Can accommodate either 500 ml or 1200 ml delivery sets
Dimensions	17" (43.2 cm) H x 8" (20.3 cm) W x 4" (10.2 cm) D

## Enteral Accessories

### Pole Clamp

FOR BOTH ENTERALITE INFINITY AND ENTERALITE ENTERAL FEEDING PUMPS




Multi-position pole clamp allows operation of pump in any orientation on poles, bed rails, and wheelchairs.

#### Specifications

Order Numbers	Z-11981
Feeding Pump Compatibility	Works with either EnteraLite Infinity or EnteraLite enteral feeding pumps

### Adjustable Angle Pole Clamp

FOR BOTH ENTERALITE INFINITY AND ENTERALITE ENTERAL FEEDING PUMPS



Multi-position and pole clamp allows operation of pump in any orientation and at nearly any angle on poles, bed rails, and wheelchairs.

#### Specifications

Order Numbers	26772-001
Feeding Pump Compatibility	Works with either EnteraLite Infinity or EnteraLite enteral feeding pumps

## Enteral Accessories

### AC Adapter/Charger with Power Cord

FOR ENTERALITE INFINITY ENTERAL FEEDING PUMP



Compact wall charger recharges pump in only 5 hours for full 24 hours of use.

#### Specifications

Order Numbers	23401-001
Feeding Pump Compatibility	Works with EnteralLite Infinity enteral feeding pump only